



## 6-Acetylmorphine Standard solution (DWT99)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	100 ug/mL and 1.0 mg/mL analytical 6-Acetylmorphine Standard
-------------------------	--

<b>Conjugate</b>	N/A
------------------	-----

<b>Cas.No</b>	2784-73-8
---------------	-----------

<b>Molecular Weight</b>	327.37g/mol
-------------------------	-------------

<b>Format</b>	Liquid
---------------	--------

<b>Concentration</b>	100 ug/mL; 1.0 mg/mL
----------------------	----------------------

<b>Size</b>	1 ml
-------------	------

<b>Buffer</b>	Acetonitrile
---------------	--------------

<b>Storage</b>	Store unopened in freezer(-10 °C to 25 °C)
----------------	--

### BACKGROUND

<b>Introduction</b>	6-Monoacetylmorphine is one of three active metabolites of heroin (diacetylmorphine), the others being morphine and the much less active 3-monoacetylmorphine (3-MAM).
---------------------	--

<b>Keywords</b>	6-MAM; 6-acetylmorphine; 6-AM; 6-Monoacetylmorphine
-----------------	---