



User's Manual

Human anti-TP ELISA Kit



IVDEIA003



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human anti-TP ELISA Kit is an enzyme-linked immunosorbent assay for qualitative determination of antibodies to *Treponema Pallidum* in human serum or plasma samples. The assay is intended to be used in clinical laboratories for diagnosis and management of patients related to infection with *Treponema Pallidum*.

Storage

The components of the kit will remain stable through the expiration date indicated on the label and package when stored at 2-8°C, do not freeze. For more detailed information, please download the following document on our website.