



# Human Anti-human GHb/HbA1c monoclonal antibody, clone K02 (CABT-L6226)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Human Anti-human GHb/HbA1c monoclonal antibody for ELISA
<b>Target</b>	Human GHb/HbA1c
<b>Immunogen</b>	Recombinant Protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	K02
<b>Purification</b>	Protein A purified. Purity ≥95% as determined by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid, Purified
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS or Tris-Gly

<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
<b>Ship</b>	Wet ice