



# HPV 66 gDNA Standard Reference (NA-GWP-SR-13)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Storage</b>	Store at 2-8°C.
<b>Product Type</b>	Nucleic Acid
<b>Buffer</b>	Tris-EDTA(10 mM Tris-HCl, 1 mM EDTA), pH 8.0
<b>Shelf Life</b>	36 months
<b>Application Note</b>	It is suitable for HPV genotyping and nucleic acid detection kit performance verification, production quality control. It is suitable for indoor quality control of HPV genotyping and nucleic acid testing, monitoring the precision and accuracy of the detection process.
<b>Concentration</b>	$2 \times 10^7$ copies/mL
<b>Unit Size</b>	100 $\mu$ L
<b>Virus type</b>	HPV 66