



HPV 84 gDNA Standard Reference (NA-GWP-SR-37)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Storage	Store at 2-8°C.
Product Type	Nucleic Acid
Buffer	Tris-EDTA(10 mM Tris-HCl, 1 mM EDTA), pH 8.0
Shelf Life	36 months
Application Note	It is suitable for HPV genotyping and nucleic acid detection kit performance verification, production quality control. It is suitable for indoor quality control of HPV genotyping and nucleic acid testing, monitoring the precision and accuracy of the detection process.
Concentration	1*10 ⁷ copies/mL
Unit Size	100 µL
Virus type	HPV 84