



Cas14a1 Nuclease (mCCE-TX-6)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Storage	Seal and store at -20°C.
Product Type	Raw Material Enzyme
Buffer	50% glycerol.
Purity	~95%(verified by SDS-PAGE)
Application Note	It can be used in multiple research areas, such as viruses and diagnosis of microbial infection, and cancer detection and modification.
Unit Size	100 pmol, 1000 pmol
Molecular Weight	61.5 kDa