



Anti-ZNF253 monoclonal antibody, clone 2E8 (CABT-14190MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description

Mouse monoclonal antibody raised against a partial recombinant ZNF253.

Immunogen ZNF253 (NP_066385, 91 a.a. ~ 180 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Isotype IgG2a

Source/Host Mouse

Species Reactivity Human

Clone 2E8

Conjugate Unconjugated

Applications sELISA, ELISA

Sequence Similarities QNSFQIGMLRRYEECRHDNLQLKKGCKSVGEHKVHKGGYNGLNQCLTTTQKEIFQCDKYG
KVFHKFSNSNTYKTRHTGINLFKCIICGKA*

Size 100 µg

Buffer In 1x PBS, pH 7.2

Preservative None

Storage Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ZNF253 zinc finger protein 253 [Homo sapiens]
Official Symbol	ZNF253
Synonyms	zinc finger protein 253; Bone marrow zinc finger 1; Zinc finger protein 411; BMZF1; BMZF-1; ZNF411; FLJ90391; DNA-binding protein
Entrez Gene ID	56242
Protein Refseq	NP_066385
UniProt ID	O75346
Chromosome Location	19p13.11
Pathway	Gene Expression, organism-specific biosystem; Generic Transcription Pathway, organism-specific biosystem.
Function	DNA binding; metal ion binding; zinc ion binding