



Anti-Rat IgD monoclonal antibody, clone MARD-3 [Biotin] (CABT-49900MR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse anti Rat IgD antibody, clone MARD-3 recognizes rat IgD, and does not cross-react with other rat immunoglobulin classes.
Specificity	IGD
Immunogen	IR731 myeloma protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Rat
Clone	MARD-3
Conjugate	Biotin
Applications	ELISA
Format	Purified IgG conjugated to Biotin - liquid
Size	500 µg
Preservative	See individual product datasheet
Storage	in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Igd immunoglobulin delta heavy chain constant region [Rattus norvegicus (Norway rat)]
Official Symbol	Igd
Synonyms	Igd; immunoglobulin delta heavy chain constant region; Ig delta chain C region
Entrez Gene ID	641523
UniProt ID	P01883
Chromosome Location	6