



Hi-Affi™ Anti-G. graminis polyclonal antibody (CABT-B1030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Gaeumannomyces graminis
Immunogen	Gaeumannomyces graminis var. graminis
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	G. graminis
Purification	Protein A affinity purified
Conjugate	Unconjugated
Applications	ELISA
Format	Lyophilized
Size	200 µg, 1 mg
Preservative	None
Storage	4°C after reconstitution.

BACKGROUND

Introduction	Gaeumannomyces graminis is a pathogenic fungi for wheat, barley, rye, oats and turf grass. Since the pathogen takes all of the grain and leaves an obvious 'white head' that is empty of seed, the disease it produces is commonly known as Take-all.
--------------	---

Keywords G. graminis;Gaeumannomyces graminis

GENE INFORMATION

Synonyms G. graminis;Gaeumannomyces graminis
