



Magic™ Anti-Asian Sea bass IgM monoclonal antibody [HRP] (CABT-B1062)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This monoclonal antibody (Mab) reacts with Asian Sea bass (<i>Lates calcarifer</i>) immunoglobulin M (IgM) and recognises the heavy chain of the IgM molecule.
Target	Asian Sea bass (<i>Lates calcarifer</i>) IgM
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Asian Sea bass
Conjugate	HRP
Applications	ELISA
Format	Lyophilized
Buffer	PBS
Preservative	None
Storage	Store at -20°C prior to reconstitution. For prolonged storage, the Mab solution should be stored at -20°C as working aliquots. Repeated freeze/thawing of the product should be avoided.

BACKGROUND

Keywords Asian Sea bass; *Lates calcarifer*

GENE INFORMATION

Synonyms	Asian Sea bass;Lates calcarifer
----------	---------------------------------
