



# Anti-CHIKV 181/25 polyclonal Antibody (CABT-B1086)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Previous lot of the antiserum showed some cross-reactivity with CHIKV OPY-1, CHIKV S27, fetal bovine serum, Vero cells, and Venezuelan Equine Encephalitis strain TC-83.
<b>Target</b>	Chikungunya Virus
<b>Immunogen</b>	181/25 virus strain
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	CHIKV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Concentration</b>	Liquid
<b>Size</b>	500 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	It is recommended that aliquots be made and stored at -20°C long term; if storing at 4°C is desired longer term, it is recommended that a preservative be added.