



Anti-CHIKV E2 polyclonal Antibody (CABT-B1087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reactivity is observed with Venezuelan Equine Encephalitis strain TC-83, Japanese Encephalitis Virus 14-14-2, or Junin Candid #1.
Target	Chikungunya Virus
Immunogen	A peptide sequence specific to Chikungunya E2 glycoprotein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	CHIKV
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	2.56 mg/mL
Size	500 µg
Buffer	PBS
Preservative	None
Storage	It is recommended that aliquots be made and stored at -20°C long term; if storing at 4°C is desired longer term, it is recommended that a preservative be added.