



Mouse anti-Human PLCL2 monoclonal antibody, clone 3E21 (CABT-B11007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PLCL2 (AAH36392, 121 a.a. ~ 211 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3E21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	RYLISYGKHTLDMLESSQDNMRTSWVSQMFSEIDVDNLGHITLCNAVQCIRNLNPGLKTS KIELKFKELEHKSKDKAGTEVTKEEFIEVFH*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	PLCL2; phospholipase C-like 2; PLCE2; inactive phospholipase C-like protein 2; PLC-L2; PLC-L(2); PLC-epsilon-2; phospholipase C-L2; phospholipase C-epsilon-2; phospholipase C, epsilon
-----------------	---

GENE INFORMATION

Entrez Gene ID	23228
----------------	-----------------------

UniProt ID	Q9UPR0
------------	------------------------

Function	calcium ion binding; phosphatidylinositol phospholipase C activity; signal transducer activity
----------	--
