



Mouse anti-Human PLEKHA6 monoclonal antibody, clone 3D4 (CABT-B11009)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PLEKHA6 (AAH10522, 1 a.a. ~ 31 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3D4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MHPRWAARLPLFISLLERADSVTAAYAKQH*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Pleckstrin homology domain-containing family A member 6 is a protein that in humans is encoded by the PLEKHA6 gene. Mouse monoclonal antibody raised against a full length recombinant PLEKHA6.
---------------------	---

Keywords

PLEKHA6; pleckstrin homology domain containing, family A member 6; PEPP3; PEPP-3; pleckstrin homology domain-containing family A member 6; PH domain-containing family A member 6; phosphoinositol 3-phosphate-binding protein 3; phosphoinositol 3-phosphate-binding protein-3;

GENE INFORMATION

Entrez Gene ID [22874](#)

UniProt ID [Q9Y2H5](#)
