



Human anti-Human TREML2 recombinant monoclonal antibody, clone SFB625 [PE] (CABT-B12152)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	SFB625
Conjugate	PE
Applications	FC
Format	Liquid
Size	100 tests
Buffer	Reagents are supplied in buffer containing stabilizer and 0.05% sodium azide.
Storage	Store protected from light at 2–8 °C. Do not freeze.

BACKGROUND

Introduction	It recognizes the human Trem-like transcript 2 (TLT-2) antigen, a single-pass type I membrane protein, which is also known as triggering receptor expressed on myeloid cells-like protein 2 (TREML2). TLT-2 is one of the conserved TREM family members, which is expressed on B cells, T cells, macrophages, and neutrophils. It is suggested to play an important role in the immune system, for it is not only highly expressed in inflammation, but also promotes the
--------------	---

proliferation of T cells and the production of cytokines. Additional information: It displays negligible binding to Fc receptors

Keywords

TREML2; triggering receptor expressed on myeloid cells-like 2; TLT2; TLT-2; C6orf76; dJ238O23.1; trem-like transcript 2 protein; triggering receptor expressed on myeloid cells-like protein 2;

GENE INFORMATION

Entrez Gene ID

[79865](#)

UniProt ID

[Q5T2D2](#)
