

Anti-DPP4 (Extracellular) monoclonal antibody, clone 0G2.3 (CABT-B1429)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with both membrane-bound (MDPP) and soluble (SDPP) forms.
Immunogen	GST-tagged recombinant protein corresponding to the extracellular domain of human DDPIV/CD26.
Isotype	lgG2a, к
Source/Host	Mouse
Species Reactivity	Human
Clone	0G2.3
Purification	Protein G Purified
Conjugate	Unconjugated
Applications	WB
Epitope	extracellular domain
Molecular Weight	~100 kDa observed. Due to glycosylation, target band appears larger than the calculated molecular weight of 88.3 (MDPP) & 84.4 (SDPP) kDa. Uncharacterized band(s) may appear in some lysates.
Format	Liquid
Concentration	Please refer to lot specific datasheet.
Size	100 µg

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Email: info@creative-diagnostics.com

Buffer	0.1 M Tris-Glycine (pH 7.4), 150 mM NaCl with 0.05% sodium azide.
Preservative	0.05% Sodium Azide
Storage	Stable for 1 year at 2-8°C from date of receipt.

BACKGROUND

Introduction

Dipeptidyl peptidase 4 (EC 3.4.14.5; UniProt P27487; also known as ADABP, ADCP-2, Adenosine deaminase complexing protein 2, DPP IV, T-cell activation antigen CD26, Adenosine deaminase complexing protein 2, Dipeptidyl peptidase 4, Dipeptidylpeptidase IV, TP103) is encoded by the DPP4 (also known as CD26, ADCP2) gene (Gene ID 1803) in human. DPPIV is a homodimeric serine exopeptidase that cleaves x-proline dipeptides from the N-terminus of polypeptides. It is involved in many cellular processes such as activation of cytokines, differentiation, and cell-matrix interactions. Inhibition of DPPIV has been reported to be an effective treatment for type II diabetes. Also, DPPIV acts as a functional receptor for the newly discovered Middle East respiratory syndrome (MERS) coronavirus.

GENE INFORMATION

Entrez Gene ID	<u>1803</u>
UniProt ID	<u>P27487</u>

© Creative Diagnostics All Rights Reserved