



Magic™ Anti-Botulinum Neurotoxin, Type A polyclonal antibody (CABT-B1972)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Botulinum Neurotoxin, Type A
Immunogen	Heavy chain domain of Botulinum Neurotoxin Type A
Isotype	IgY
Source/Host	Chicken
Species Reactivity	C. botulinum
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, sELISA
Concentration	Store this product at 2-8°C
Size	100 µg
Buffer	PBS, pH 7.4 with 0.075% azide
Preservative	0.075% Sodium Azide
Storage	Store this product at 2-8°C