



Magic™ Anti-Botulinum Neurotoxin, Type A monoclonal antibody (CABT-B1973)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Target | Botulinum Neurotoxin, Type A |
| Immunogen | Light chain of Botulinum Neurotoxin type A |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | C. botulinum |
| Conjugate | Unconjugated |
| Applications | WB, sELISA, HPLC |
| Concentration | Store this product at 2-8°C. Once dissolved, aliquot and store the product at -20°C. Refrain from multiple freeze/thaw cycles. |
| Size | 100 µg |
| Buffer | PBS, pH 7.4 with 0.05% azide |
| Preservative | 0.05% Sodium Azide |
| Storage | Store this product at 2-8°C. Once dissolved, aliquot and store the product at -20°C. Refrain from multiple freeze/thaw cycles. |