



Anti-P. aeruginosa Exotoxin A polyclonal antibody (CABT-B1977)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	P. aeruginosa Exotoxin A
Immunogen	Exotoxin A from Pseudomonas aeruginosa
Source/Host	Goat
Species Reactivity	P. aeruginosa
Conjugate	Unconjugated
Applications	ELISA, WB
Reconstitution	0.5 ml water
Concentration	Store lyophilized vials at 2-8°C. Once dissolved, aliquot and store the product at -20°C. Refrain from multiple freeze/thaw cycles.
Size	0.5 ml
Preservative	None
Storage	Store lyophilized vials at 2-8°C. Once dissolved, aliquot and store the product at -20°C. Refrain from multiple freeze/thaw cycles.