



Anti-Llama IgG (H&L) polyclonal antibody (CABT-B1991)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Based on IEP, this antibody reacts with: heavy (γ) chains on llama IgG light chains on all llama immunoglobulins
Immunogen	Purified Llama IgG, whole molecule
Source/Host	Goat
Species Reactivity	Llama
Purification	Affinity purified using solid phase Llama IgG
Conjugate	Unconjugated
Applications	IA
Reconstitution	Rehydrate with 1.1 ml of deionized water and let stand 30 minutes at room temperature to dissolve. Centrifuge to remove any particulates. Prepare fresh working dilution daily.
Format	Clear, colorless liquid, 0.2 um filtered
Concentration	> 4.5 mg/ml
Size	1 mg
Buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative	0.05% Sodium Azide
Storage	2-8°C