



Anti-Llama IgG Light Chain polyclonal antibody (CABT-B2010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified Llama IgG Light Chain
Source/Host	Rabbit
Species Reactivity	Llama
Purification	Affinity Purification
Conjugate	Unconjugated
Applications	WB
Size	1 mg
Buffer	PBS (Phosphate Buffered Saline) pH: 7.4 with 0.02% azide
Preservative	0.02% Sodium Azide
Storage	Aliquot and freeze at -20°C. For long term storage freeze at -70°C. Avoid freeze-thaw cycles.