



Anti-Rabbit IgG polyclonal antibody (CABT-B2019)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified Rabbit IgG
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Rabbit
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid, purified
Size	1 mg
Buffer	Tris buffer with 0.1% sodium azide as a preservative.
Preservative	0.1% Sodium Azide
Storage	Aliquot and freeze at -20°C. Avoid freeze-thaw cycles.