



Anti-Chicken IgY polyclonal antibody [FITC] (CABT-B2021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	purified Chicken IgY
Source/Host	Rabbit
Species Reactivity	Chicken
Conjugate	FITC
Applications	FC
Format	Liquid, purified
Size	100 μΙ
Buffer	Phosphate buffer, pH 7.4
Preservative	None
Storage	Aliquot and freeze at -20°C. Avoid freeze-thaw cycles.

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221