



Anti-ACBD3 (Internal) polyclonal antibody (CABT-B361)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	ACBD3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Dog, Rat, Zebrafish, Yeast, Guinea pig, Horse, Rabbit, Bovine
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration)
Size	100 µl
Buffer	In 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Preservative	0.09% Sodium Azide
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.

BACKGROUND

Keywords ACBD3;acyl-CoA binding domain containing 3;PAP7;GCP60;GOCAP1;GOLPH1;Golgi resident

protein GCP60;golgi phosphoprotein 1;PKA (RIalpha)-associated protein;PBR- and PKA-associated protein 7;Acbd3

GENE INFORMATION

Synonyms acyl-CoA binding domain containing 3; Acbd3
