



Anti-AKT3 polyclonal antibody (CABT-B384)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	AKT3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Purification	Purified by antigen-affinity chromatography.
Conjugate	Unconjugated
Applications	IHC-Wmt
Positive Control	Zebrafish embryo
Format	Liquid
Concentration	0.78 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µl
Buffer	1XPBS, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Preservative	None
Storage	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Keywords AKT3;v-akt murine thymoma viral oncogene homolog 3 (protein kinase B, gamma);RAC-

gamma serine/threonine-protein kinase;PKBG;PRKBG;RAC gamma;PKB gamma;RAC-gamma
serine/threonine protein kinase;STK-2;PKB-GAMMA;RAC-gamma;RAC-PK-
gamma;DKFZp434N0250

GENE INFORMATION

Synonyms RAC-gamma serine/threonine-protein kinase; v-akt murine thymoma viral oncogene homolog 3
(protein kinase B, gamma), akt3
