



Anti-ASCL2 polyclonal antibody (CABT-B401)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Target | ASCL2 |
| Immunogen | Synthetic peptide corresponding to a region of Mouse |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human, Mouse, Dog, Rat, Zebrafish, Guinea pig, Rabbit, Bovine |
| Conjugate | Unconjugated |
| Applications | WB |
| Format | Lyophilized powder |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration) |
| Size | 100 µg |
| Buffer | PBS buffer with 2% sucrose. |
| Preservative | None |
| Storage | For longer periods of storage, store at -20°C. Avoid repeat freeze-thaw cycles. |

BACKGROUND

| | |
|-----------------|--|
| Keywords | ASCL2;achaete-scute complex homolog 2 (Drosophila);achaete scute complex (Drosophila) homolog like 2;achaete scute complex like 2 (Drosophila);achaete-scute homolog 2;ASH2;bHLHa45;HASH2;ASH-2;achaete-scute complex-like 2 |
|-----------------|--|

GENE INFORMATION

Synonyms

bHLHa45; Mash2; achaete-scute complex homolog 2 (Drosophila); Ascl2; 2410083I15Rik
