



Anti-CDX4 (Internal) polyclonal antibody (CABT-B464)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CDX4
Immunogen	Synthetic peptide directed towards the middle region of human CDX4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Dog, Pig, Rat, Zebrafish, Guinea pig, Horse, Rabbit, Bovine
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilized powder
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µg
Buffer	PBS buffer with 2% sucrose.
Preservative	None
Storage	For longer periods of storage, store at -20°C. Avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction	This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior
---------------------	--

patterning and hematopoiesis.

Keywords

CDX4;caudal type homeobox 4;homeobox protein CDX-4;caudal-type homeobox protein 4;caudal type homeobox transcription factor 4

GENE INFORMATION

Synonyms

caudal type homeobox 4
