



Anti-Commd2 polyclonal antibody (CABT-B486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Commd2
Immunogen	Synthetic peptide corresponding to a region of Rat
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Dog, Pig, Rat, Zebrafish, Guinea pig, Horse, Rabbit, Bovine
Conjugate	Unconjugated
Applications	WB
Format	Lyophilized powder
Concentration	Lot specific
Size	50 µg
Buffer	PBS buffer with 2% sucrose.
Preservative	None
Storage	For longer periods of storage, store at -20°C. Avoid repeat freeze-thaw cycles.

BACKGROUND

Keywords	COMMD2;COMM domain containing 2;COMM domain-containing protein 2;HSPC042;MGC57611
----------	---

GENE INFORMATION

Synonyms	COMM domain containing 2
----------	--------------------------
