



## Anti-Nap1l1 polyclonal antibody (CABT-B755)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	Nap1l1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish, Human, Mouse, Rat, Bovine
<b>Purification</b>	Purified by antigen-affinity chromatography.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Positive Control</b>	Zebrafish eye
<b>Format</b>	Liquid
<b>Concentration</b>	0.92 mg/ml (Please refer to the vial label for the specific concentration)
<b>Size</b>	50 µl
<b>Buffer</b>	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
<b>Preservative</b>	None
<b>Storage</b>	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

### BACKGROUND

**Keywords** NAP1L1;nucleosome assembly protein 1-like 1;NRP;NAP1;NAP1L;hNRP;NAP-1 related

protein;NAP-1-related protein;HSP22-like protein interacting protein;

---

## GENE INFORMATION

### Synonyms

nap111b; nucleosome assembly protein 1-like 1; cb716; sb:cb13; wu:fb83g07; nap111;  
nucleosome assembly protein 1, like 1

---