



Anti-Nfyba polyclonal antibody (CABT-B760)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Nfyba
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Purification	Purified by antigen-affinity chromatography.
Conjugate	Unconjugated
Applications	IHC-P, WB, IHC-Wmt
Positive Control	Zebrafish eye, Zebrafish embryo, Ubiquitous express
Format	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µl
Buffer	1XPBS, 1%BSA, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Preservative	None
Storage	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Keywords NFYB;nuclear transcription factor Y, beta;nuclear transcription factor Y subunit beta;CBF

A;HAP3;NF YB;Transcription factor NF-Y, B subunit;CAAT box DNA-binding protein subunit B;nuclear transcription factor Y subunit B;CCAAT-binding transcription factor subunit A;CBF-A;CBF-B;NF-YB;

GENE INFORMATION

Synonyms nuclear transcription factor Y, beta; zgc:110552; nfyba; nuclear transcription factor Y, beta a
