



Anti-Ptp4a3 polyclonal antibody (CABT-B816)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Ptp4a3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish, Human, Mouse, Monkey, Rat, Bovine
Purification	Purified by antigen-affinity chromatography.
Conjugate	Unconjugated
Applications	WB
Positive Control	Zebrafish eye
Format	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µl
Buffer	1XPBS, 1%BSA, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Preservative	None
Storage	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Keywords fc54b05;fv52d11;protein tyrosine phosphatase type IVA

3;wu:fc54b05;wu:fv52d11;zgc:77109;ptp4a3;protein tyrosine phosphatase type IVA, member 3;protein tyrosine phosphatase type IVA, member 3

GENE INFORMATION

Synonyms fc54b05; fv52d11; protein tyrosine phosphatase type IVA 3; wu:fc54b05; wu:fv52d11; zgc:77109; ptp4a3; protein tyrosine phosphatase type IVA, member 3
