



Anti-Racgap1 polyclonal antibody (CABT-B821)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Racgap1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Purification	Purified by antigen-affinity chromatography.
Conjugate	Unconjugated
Applications	IHC-Wmt
Positive Control	Zebrafish embryo
Format	Liquid
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration)
Size	50 µl
Buffer	1XPBS, 1%BSA, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Preservative	None
Storage	Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Keywords RACGAP1;Rac GTPase activating protein 1;rac GTPase-activating protein

1;MgcRacGAP;protein CYK4 homolg;male germ cell RacGap;GTPase activating protein;CYK4;ID-GAP;HsCYK-4;wu:fb92c08;zgc:77839;racgap1

GENE INFORMATION

Synonyms rac GTPase-activating protein 1; wu:fb92c08; zgc:77839; racgap1
