



Anti-Zika Virus NS1 monoclonal antibody (CABT-B8492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with Zika Virus NS1
Immunogen	Zika Virus NS1 Recombinant
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	zika virus
Purification	DEAE Chromatography
Conjugate	Unconjugated
Applications	ELISA(cap), LFIA
Format	Liquid
Concentration	10.00 mg/mL
Size	1 mg
Buffer	0.015 M Potassium Phosphate Buffer, pH 7.2, 0.85% Sodium Chloride.
Preservative	0.05% Sodium Azide
Storage	Short term store at 2-8°C. For long term storage aliquot and store at -20°C to avoid repeated freeze-thaw cycles.

BACKGROUND

Keywords

Zika virus; ZIKV; Flaviviridae; Flavivirus; ZIKV NS1; Zika virus NS1; Zikv nonstructural protein 1
