



Anti-Zika Virus NS1 polyclonal antibody (CABT-B8549)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Antibody to Zika Virus Nonstructural Protein 1 (NS1), Affinity Purified
Specificity	Reacts with Zika Virus NS1.
Immunogen	Full-length Zika Virus NS1 Recombinant
Source/Host	Goat
Species Reactivity	Zika virus
Purification	> 85% pure (SDS-PAGE). Immunoaffinity purified by Zika Virus NS1 conjugated Sepharose 4B matrix
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C