



Anti-CHIKV Env monoclonal antibody, clone IN449 (CABT-B8658)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No reactivity with Dengue and Zika viruses
Immunogen	Recombinant Chikungunya Virus Envelope Protein
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	CHIKV
Clone	IN449
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Preservative	0.05% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.