



Anti-CHIKV monoclonal antibody, clone N97701 (CABT-B8659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Does not react with EEE, WEE, VEE or other Arboviruses.
Immunogen	Chikungunya virus antibody was raised in Mouse using a Chikungunya virus lysate
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	CHIKV
Clone	N97701
Purification	Protein A chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Buffer	Supplied in liquid form in 0.01 M PBS, pH 7.2 with 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 2-8°C for short term storage. Aliquot and store at -20°C for long term storage.