



Anti-NDV RNP monoclonal antibody (CABT-B8665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Ribonucleoprotein (RNP) of Newcastle Disease Virus
Immunogen	Newcastle disease virus (La-sota strain).
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	NDV
Purification	Protein G Sepharose chromatography, $\geq 95\%$ pure (SDS-PAGE).
Conjugate	Unconjugated
Applications	ELISA, IHC, IF
Format	Purified, Liquid
Concentration	7.7 mg/mL(OD280nm, E0.1% = 1.4)
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C.