



# Magic™ Anti-E. coli Enterotoxin LT monoclonal antibody, clone C2032N (CABT-B8711)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Specific to E. coli Heat-Labile Enterotoxin (LT). Cross reactive to Cholera toxin (CT).
<b>Immunogen</b>	Heat-Labile Enterotoxin
<b>Isotype</b>	IgG2a
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	E. coli
<b>Clone</b>	C2032N
<b>Purification</b>	Protein A Chromatography
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA (cap)
<b>Format</b>	Liquid
<b>Buffer</b>	PBS: 1.5 mM Monopotassium Phosphate (KH <sub>2</sub> PO <sub>4</sub> ), 8.0 mM Disodium Phosphate (Na <sub>2</sub> HPO <sub>4</sub> ), 3 mM Potassium Chloride (KCl), 140 mM Sodium Chloride (NaCl), pH 7.4.
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.