



Magic™ Anti-E. coli Enterotoxin LT monoclonal antibody, clone C2035N (CABT-B8712)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to E. coli Heat-Labile Enterotoxin (LT). Cross reactive to Cholera toxin (CT).
Immunogen	Heat-Labile Enterotoxin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	E. coli
Clone	C2035N
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA (det)
Format	Liquid
Buffer	PBS: 1.5 mM Monopotassium Phosphate (KH ₂ PO ₄), 8.0 mM Disodium Phosphate (Na ₂ HPO ₄), 3 mM Potassium Chloride (KCl), 140 mM Sodium Chloride (NaCl), pH 7.4.
Preservative	0.05% Sodium Azide
Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.