



Anti-*S. pneumoniae* common C-polysaccharide monoclonal antibody, clone C2400N (CABT-B8728)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reactive with the <i>Streptococcus pneumoniae</i> common C-polysaccharide, phosphorylcholine epitope.
Immunogen	Purified <i>Strep. pneumoniae</i> common cell wall polysaccharide.
Isotype	IgG3
Source/Host	Mouse
Species Reactivity	<i>S. pneumoniae</i>
Clone	C2400N
Purification	Protein A Chromatograph
Conjugate	Unconjugated
Applications	ELISA, IF
Format	Liquid
Size	1 mg
Buffer	0.01 M PBS, pH 7.2 Product contains no stabilizing proteins
Preservative	0.1% Sodium Azide
Storage	Upon receipt, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.