



Anti-C. parvum (Oocysts) monoclonal antibody, cloen PX61H22 (CABT-B8730)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cryptosporidium species (oocysts).
Isotype	IgG3
Source/Host	Mouse
Species Reactivity	Cryptosporidium parvum
Clone	PX61H22
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	IF, LFIA
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4 + 0.2
Preservative	0.05% Sodium Azide
Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.