



# Anti-*C. parvum* (Oocysts) monoclonal antibody, clone CEJ481 (CABT-B8731)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	<i>C. parvum</i> .
Immunogen	Purified oocysts.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	<i>Cryptosporidium parvum</i>
Clone	CEJ481
Purification	Protein A Chromatography
Conjugate	Unconjugated
Applications	IF, IHC, ELISA
Format	Liquid
Size	1 mg
Buffer	0.01 M PBS, pH 7.2 Product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Storage	Upon receipt, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles