



## Anti-HPV11 L1 polyclonal antibody [HRP] (CABT-B8788)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HPV11
Purification	Protein G purified
Conjugate	HRP
Applications	WB, ELISA
Format	Liquid
Concentration	1mg/1ml
Size	100 µg
Buffer	10mM PBS, 25% glycerol
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing.