



Anti-HPV18 L1 polyclonal antibody (CABT-B8791)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HPV18
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	1mg/1ml
Size	100 µg
Buffer	Tris-Gly
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing.