



# Anti-HPV 11 L1 monoclonal antibody (set) (CABT-B8804)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	The antibody set can be used for ELISA and screening appropriate mAbs to establish VLP vaccine in vitro specific activity assay
<b>Immunogen</b>	E. coli expressed HPV11 L1 VLP
<b>Isotype</b>	34C9-IgG1; 35C1-IgG1; 35H5-IgG1; 34E5-IgG1; 40C5-IgG1; 34D10-IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	HPV 11
<b>Clone</b>	34C9; 35C1; 35H5; 34E5; 40C5; 34D10
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, Screening
<b>Format</b>	Liquid
<b>Concentration</b>	0.5mg/ml
<b>Size</b>	20 µgx6
<b>Buffer</b>	Tris-Gly
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Avoid repeated freezing and thawing.