



Anti-HPV 45 L1 monoclonal antibody (set) (CABT-B8809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The antibody set can be used for ELISA and screening appropriate mAbs to establish VLP vaccine in vitro specific activity assay
Immunogen	E. coli expressed HPV45 L1 VLP
Isotype	30H6-IgG1; 35B12-IgG2b; 35H9-IgG1; 40H1-IgG2a; 46G5-IgG2a; 48B1-IgG2a
Source/Host	Mouse
Species Reactivity	HPV 45
Clone	30H6; 35B12; 35H9; 40H1; 46G5; 48B1
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, Screening
Format	Liquid
Concentration	0.5mg/ml
Size	20 µgx6
Buffer	Tris-Gly
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing.