



Anti-Human IgG polyclonal antibody [FITC] (CABT-B8937)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Specificity | Human IgG specific |
| Immunogen | Whole molecule human IgG |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunoaffinity chromatography |
| Conjugate | FITC |
| Applications | IHC, FC, ICC |
| Format | Liquid |
| Size | 1 ml |
| Buffer | 0.01 M PBS (PH 7.4), 0.01% Thimerosal, 50% glycerol |
| Preservative | None |
| Storage | 4°C, store in dark |