



## Anti-DM1 polyclonal antibody (CABT-B8996)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Target	DM1
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Purification	Affinity Chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	This antibody comes as 50x concentrated and is used in the DM1-ADC competitive EIA
Size	100 µl
Buffer	0.01M PBS, 0.1%BSA, Proclin-300
Preservative	None
Storage	-20°C (-80°C for long term)

### BACKGROUND

Keywords	DM1
----------	-----