



## RHA™ Anti-Trenbolone monoclonal antibody, clone TB (CABT-B9092)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	Trenbolone-BSA
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TB
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Shipped at 4°C

### BACKGROUND

**Introduction** Trenbolone (INN, BAN), also known as trienolone or trienbolone, as well as 11-testosterone, is an anabolic-androgenic steroid (AAS) of the 19-nortestosterone group. It is used in veterinary medicine in livestock to increase muscle growth and appetite. Trenbolone is used in the form of ester prodrugs, including trenbolone acetate (brand names Revalor (in combination with estradiol), Finaplix, Finajet) and trenbolone hexahydrobenzylcarbonate (Parabolan, Hexabolan).Although it is not available as a prescription medication for humans, trenbolone enanthate (RU-1697) is sometimes sold under the nickname Trenabol

<b>Keywords</b>	Trenbolone
-----------------	------------

## GENE INFORMATION

<b>Gene Name</b>	Trenbolone
<b>Synonyms</b>	Trenbolone