



Mouse anti-Human ARL9 monoclonal antibody, clone 8B3 (CABT-B9793)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ARL9 (NP_996802.1, 1 a.a. ~ 123 a.a) full-length recombinant protein with GST tag.
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Human
Clone	8B3
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	MEFLEIGGSKPFRSYWEMYLSKGLLLIFVVDSADHSRLPEAKKYLHQLIAANPVLPLVVF ANKQDLEAAYHITDIHEALALSEVGNDRKMFLFGTYLTKNGSEIPSTMQDAKDLIAQLAA DVQ
Format	Liquid
Size	200 µl
Buffer	In ascites fluid
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	ARL9 is a member of the small GTPase protein family with a high degree of similarity to ARF (MIM 103180) proteins of the RAS superfamily.[supplied by OMIM, Nov 2008]
---------------------	---

Keywords

ARL9; ADP-ribosylation factor-like 9; ADP-ribosylation factor-like protein 9;

GENE INFORMATION

Entrez Gene ID[132946](#)

UniProt ID[Q6T311](#)

FunctionGTP binding; nucleotide binding
